

UniDyn--Demo-07.nb

Aritro Sinha Roy
as836@cornell.edu
Cornell University

Abstract: This demonstration notebook loads the **UniDyn** package and illustrate its use in simulating a three pulse experiment with phase-cycling for selective generation of specific quantum coherence orders using a system of three spin-1/2 particles.

Set the path to the package

Tell *Mathematica* the path to the directory containing the package.

EDIT THE FOLLOWING PATH STRING:

```
In[43]:= $UniDynPath =  
        "C:/Users/as836/Documents/Work_with_John/UniDyn-master/unidyn";
```

YOU SHOULD NOT NEED TO EDIT ANYTHING FROM HERE ONWARDS.

Load the package

Append the package path to the system path. Before trying to load the package, ask *Mathematica* to find it. This is a test that we directed *Mathematica* to the correct directory. The output of this command should be the full system path to the UniDyn.m file.

```
In[44]:= (* $Path = AppendTo[$Path,$NCPATH];*)  
$Path = AppendTo[$Path, $UniDynPath];  
FindFile["UniDyn`"]
```

```
Out[45]=  
C:\Users\as836\Documents\Work_with_John\UniDyn-master\unidyn\UniDyn.m
```

Now that we are confident that the path is set correctly, load the package. Setting the global \$VerboseLoad variable to True will print out the help strings for key commands in the package.

```
In[46]:= $VerboseLoad = True; (* Set to load quietly *)
Needs["UniDyn`"]
```

Run the unit tests .

```
In[48]:= SetDirectory[$UniDynPath];
test$report = TestReport /@ FileNames["*-tests.m"];
TableForm[Table[test$report[[k]], {k, 1, Length[test$report]}]]
```

Out[50]//TableForm=

TestReportObject	 	Title: Test Report: Comm-tests.m Success rate: 100% Tests count: 23
TestReportObject	 	Title: Test Report: Evolver1-tests.m Success rate: 100% Tests count: 14
TestReportObject	 	Title: Test Report: Evolver2-tests.m Success rate: 100% Tests count: 13
TestReportObject	 	Title: Test Report: Evolve-tests.m Success rate: 100% Tests count: 9
TestReportObject	 	Title: Test Report: Inv-tests.m Success rate: 100% Tests count: 24
TestReportObject	 	Title: Test Report: Mult-tests.m Success rate: 100% Tests count: 18
TestReportObject	 	Title: Test Report: OpQ-tests.m Success rate: 100% Tests count: 21
TestReportObject	 	Title: Test Report: Osc-tests.m Success rate: 100% Tests count: 22
TestReportObject	 	Title: Test Report: SpinBoson-tests.m Success rate: 100% Tests count: 8
TestReportObject	 	Title: Test Report: Spins-tests.m Success rate: 100% Tests count: 14

Three spin-1/2 particles: Selective quantum coherence generation

Create a system of three spin-1/2 particles.

```
In[51]:= Needs["UniDyn`"];
Clear[Ix1, Iy1, Iz1, Δ1, Ix2, Iy2,
      Iz2, Δ2, Ix3, Iy3, Iz3, Δ3, a12, a13, a23]
Quiet[
CreateScalar[{Δ1, Δ2, Δ3, a12, a13, a23}];
CreateOperator[{{Ix1, Iy1, Iz1}, {Ix2, Iy2, Iz2}, {Ix3, Iy3, Iz3}}];
SpinSingle$CreateOperators[Ix1, Iy1, Iz1, 1/2];
SpinSingle$CreateOperators[Ix2, Iy2, Iz2, 1/2];
SpinSingle$CreateOperators[Ix3, Iy3, Iz3, 1/2];
]
```

Define some useful functions to extract operators

```
In[54]:= (* Extraction of operators from an expression *)
op$extA[elem_] := Module[{pos, ans},
  pos = Position[OperatorQ /@ Level[elem, 1], True];
  ans = Level[elem, 1][[pos[[1, 1]]];
  Return[ans];
];

op$ext[elem_] := Module[{ans},
  If[Dimensions[elem] == {}, ans = elem,
    ans = Check[op$extA[elem], If[OperatorQ@elem, elem, None]] // Quiet];
  Return[ans];
];

(* Group and simplify expressions by operators *)
col$op[exp_, op_, s_] := Module[{op$l, exp$l, dim, exp$s, j}, exp$l = exp // Expand;
  dim = Dimensions[exp$l][[1]];
  op$l = {};
  For[j = 1, j ≤ dim, j++, op$l = Append[op$l, op$ext[exp$l][[j]]]];
  op$l = DeleteDuplicates[op$l];
  If[op == 1, exp$s = Simplify[#, TimeConstraint → s] & /@ Collect[exp$l, op$l],
    If[op == 2, exp$s = FullSimplify[#, TimeConstraint → s] & /@ Collect[exp$l, op$l],
      exp$s = Collect[exp$l, op$l]];
  Return[exp$s] ;
```

Define free-evolution and pulse propagators

```

In[57]:= (* Define free-evolution under the influence of the secular spin Hamiltonian *)
FreeEvolution[ρ_, t_, sim_] := (*ρ=density operator*)
(*t=time[s]*) (* sim = 1: Simplify, 2: FullSimplify *)
Module[{dim, A, B, Anew, ρnew},
  dim = Dimensions[ρ][[1]];
  ρnew =
    (Evolve[Δ1 Iz1, t, Evolve[Δ2 Iz2, t, Evolve[Δ3 Iz3, t, Evolve[a12 Mult[Iz1, Iz2], t,
      Evolve[a13 Mult[Iz1, Iz3], t, Evolve[a23 Mult[Iz2, Iz3], t, ρ] /.
        Evolve → Evolver2 /. Mult → SortedMult] /. Evolve →
          Evolver2 /. Mult → SortedMult] /. Evolve → Evolver2 /.
            Mult → SortedMult] /. Evolve → Evolver2 /. Mult → SortedMult] /.
              Evolve → Evolver2 /. Mult → SortedMult] /. Evolve → Evolver2 /.
                Mult → SortedMult) // Expand // MultSort // Expand;
  ρnew = col$op[ρnew, sim, 1]; (* Timeconstraint = 1 for simplification *)
  Return[ρnew];
];

(* Define non-selective ideal δ-pulses with flip angle θ and phase φ *)
Pulse[ρ_, θ_, φ_, sim_] := (*ρ=density operator*) (*θ=flip angle*) (*φ=phase*)
Module[{dim, A, B, Anew, ρnew},
  ρnew =
    (Evolve[Ix1 Cos[φ] + Iy1 Sin[φ], θ, Evolve[Ix2 Cos[φ] + Iy2 Sin[φ], θ, Evolve[Ix3 Cos[
      φ] + Iy3 Sin[φ], θ, ρ] /. Evolve → Evolver2 /. Mult → SortedMult] /.
      Evolve → Evolver2 /. Mult → SortedMult] /. Evolve → Evolver2 /.
        Mult → SortedMult) // Expand // MultSort // Expand // Quiet;
  ρnew = col$op[ρnew, sim, 1]; (* Timeconstraint = 1 for simplification *)
  Return[ρnew];
];

```


Overview of the pulse-sequence

```

In[59]:= text$1 = Text[Style["t1", Large], {2.5, 0.6}];
text$2 = Text[Style["t1", Large], {5.5, 0.6}];
gr$ = Graphics[{text$1, text$2}];
pp$ = BarChart[{1, Missing[], Missing[], 2, Missing[], Missing[], 1},
  ChartStyle → {"Pastel", LightGray}, ChartBaseStyle → EdgeForm[None],
  ChartLabels → {"π/2-pulse", "", "", "π-pulse", "", "", "π/2-pulse"},
  LabelStyle → {FontSize → 22},
  Frame → {True, True, False, False},
  FrameTicksStyle → {Opacity[0], Opacity[1]}, FrameStyle → Opacity[0],
  (* Axes → {True, False}, *) ImageSize → 400];
DeleteCases[pp$, _Line? (Not@*FreeQ[_Offset]), All];
Show[pp$, gr$]

```

Out[64]=



Density operator evolution

```

In[65]:= AbsoluteTiming[
  (* ρ0 is the density operator at thermal equilibrium *)
  ρ0 = Iz1;

  (* Application of the first π/2-pulse with its phase set to {0, π/2, π, 3 π/2} *)
  ρ1 = Pulse[ρ0, π/2, φ1, 2] // Quiet;

  (* A free-evolution period *)
  ρ2 = FreeEvolution[ρ1, t1, 2] // Quiet;

  (* Application of the first refocusing π-
  pulse with its phase set to {0, π/2, π, 3 π/2} *)
  ρ3 = 0;
  dim = Dimensions[ρ2][[1]];
  Monitor[
    For[k = 1, k ≤ dim, k++,
      ρ = Pulse[ρ2[[k]], π, φ2, 1];

```

```

AddTo[ $\rho_3$ ,  $\rho$ ]],
Row[{ProgressIndicator[k, {1, dim}], NumberForm[ $\frac{1. * k}{dim}$ , {2, 2}]}],
" %  $\rho_3$  calculated = "]];
 $\rho_{3F}$  = col$op[ $\rho_3$ , 2, 1];

(* The second free-evolution period *)
 $\rho_4$  = 0;
dim = Dimensions[ $\rho_{3F}$ ][[1]];
Monitor[
For[k = 1, k ≤ dim, k++,
 $\rho$  = FreeEvolution[ $\rho_{3F}$ [[k]], t1, 1];
AddTo[ $\rho_4$ ,  $\rho$ ]],
Row[{ProgressIndicator[k, {1, dim}], NumberForm[ $\frac{1. * k}{dim}$ , {2, 2}]}],
" %  $\rho_4$  calculated = "]];
 $\rho_{4F}$  = col$op[ $\rho_4$ , 2, 1];

(* The second  $\pi/2$ -pulse to create multi-quantum coherence *)
 $\rho_5$  = 0;
dim = Dimensions[ $\rho_{4F}$ ][[1]];
Monitor[
For[k = 1, k ≤ dim, k++,
 $\rho$  = Pulse[ $\rho_{4F}$ [[k]],  $\pi/2$ ,  $\phi_3$ , 1];
AddTo[ $\rho_5$ ,  $\rho$ ]],
Row[{ProgressIndicator[k, {1, dim}], NumberForm[ $\frac{1. * k}{dim}$ , {2, 2}]}],
" %  $\rho_5$  calculated = "]];
 $\rho_{5F}$  = col$op[ $\rho_5$ , 2, 1];][[1]]

```

Out[65]=

6.04068

The general analytical expression of the density operator in the three pulse experiment is given by ρ_{5F} . We will show how this expression can be used efficiently in understanding phase cycles that selectively generate specific coherence orders.

****Note that the initial density operator ρ_0 is taken to be I_{z1} : given that the problem is symmetric with respect to the three spins, the density operator expressions for $\rho_0 = I_{z2}$ or I_{z3} can be obtained by permutation of the indices in ρ_{5F} .****

In general, three spin-1/2 particles can generate coherence orders in the range of ± 3 . By setting $\phi_1 = \phi_2 = \phi_3$, one can selectively generate only even order coherence.

```
In[66]:=  $\rho_{\text{EvenQC}} = \text{col\$op}[\rho_{5F} /. \{\phi_2 \rightarrow \phi_1, \phi_3 \rightarrow \phi_1\}, 2, 1]$ 
```

```
Out[66]=
```

$$\begin{aligned} & I_z 1 \cos[a_{12} t_1] \cos[a_{13} t_1] + 2 \cos[a_{13} t_1] \cos[\phi_1]^2 \text{Mult}[I_{x1}, I_{y2}] \sin[a_{12} t_1] + \\ & 2 \cos[a_{12} t_1] \cos[\phi_1]^2 \text{Mult}[I_{x1}, I_{y3}] \sin[a_{13} t_1] - \\ & 4 \cos[\phi_1]^2 \text{Mult}[I_{z1}, I_{y2}, I_{y3}] \sin[a_{12} t_1] \sin[a_{13} t_1] - \\ & 2 \cos[a_{13} t_1] \text{Mult}[I_{y1}, I_{x2}] \sin[a_{12} t_1] \sin[\phi_1]^2 - \\ & 2 \cos[a_{12} t_1] \text{Mult}[I_{y1}, I_{x3}] \sin[a_{13} t_1] \sin[\phi_1]^2 - \\ & 4 \text{Mult}[I_{z1}, I_{x2}, I_{x3}] \sin[a_{12} t_1] \sin[a_{13} t_1] \sin[\phi_1]^2 - \\ & \cos[a_{13} t_1] \text{Mult}[I_{x1}, I_{x2}] \sin[a_{12} t_1] \sin[2 \phi_1] + \\ & \cos[a_{13} t_1] \text{Mult}[I_{y1}, I_{y2}] \sin[a_{12} t_1] \sin[2 \phi_1] - \\ & \cos[a_{12} t_1] \text{Mult}[I_{x1}, I_{x3}] \sin[a_{13} t_1] \sin[2 \phi_1] + \\ & \cos[a_{12} t_1] \text{Mult}[I_{y1}, I_{y3}] \sin[a_{13} t_1] \sin[2 \phi_1] + \\ & 2 \text{Mult}[I_{z1}, I_{x2}, I_{y3}] \sin[a_{12} t_1] \sin[a_{13} t_1] \sin[2 \phi_1] + \\ & 2 \text{Mult}[I_{z1}, I_{y2}, I_{x3}] \sin[a_{12} t_1] \sin[a_{13} t_1] \sin[2 \phi_1] \end{aligned}$$

One can inspect the above expression to find that it contains terms corresponding to 0 and ± 2 coherence orders only.

For example, the following terms from ρ_{EvenQC}

$$A = -4 \cos[\phi_1]^2 \text{Mult}[I_{z1}, I_{y2}, I_{y3}] \sin[a_{12} t_1] \sin[a_{13} t_1] - 4 \text{Mult}[I_{z1}, I_{x2}, I_{x3}] \sin[a_{12} t_1] \sin[a_{13} t_1] \sin[\phi_1]^2$$

and lowering operators as:

$$A = \left(\cos[\phi_1]^2 \text{Mult}[I_{z1}, (I_{p2} - I_{m2}), (I_{p3} - I_{m3})] - \text{Mult}[I_{z1}, (I_{p2} + I_{m2}), (I_{p3} + I_{m3})] \sin[\phi_1]^2 \right) \sin[a_{12} t_1] \sin[a_{13} t_1]$$

Setting $\phi_1 \rightarrow 0$ and $\phi_1 \rightarrow \pi/2$, the above expression transforms to

$$A_1 = \text{Mult}[I_{z1}, (I_{p2} - I_{m2}), (I_{p3} - I_{m3})] \sin[a_{12} t_1] \sin[a_{13} t_1], \text{ and}$$

$$A_2 = -\text{Mult}[I_{z1}, (I_{p2} + I_{m2}), (I_{p3} + I_{m3})] \sin[a_{12} t_1] \sin[a_{13} t_1]$$

Now when we add or subtract A_1 and A_2 , we can produce zero or double quantum coherence selectively, given by

$$A_1 + A_2 =$$

$$-2 \times (\text{Mult}[I_{z1}, I_{p2}, I_{m3}] + \text{Mult}[I_{z1}, I_{m2}, I_{p3}]) \sin[a_{12} t_1] \sin[a_{13} t_1]$$

$$A_1 - A_2 = 2 \times (\text{Mult}[I_{z1}, I_{p2}, I_{p3}] + \text{Mult}[I_{z1}, I_{m2}, I_{m3}]) \sin[a_{12} t_1] \sin[a_{13} t_1]$$

This is the basis of phase cycling in selecting quantum coherence of different orders in pulsed magnetic resonance experiments.

We apply the phase cycle to ρ_{EvenQC} to separate out 0 and ± 2 coherence orders

```
In[67]:= zq$ = FullSimplify[Sum[(1)^(2*j)*rhoEvenQC /. phi1 -> j * (pi/2), {j, 0, 1}]/2];
dq$ = FullSimplify[Sum[(-1)^j * rhoEvenQC /. phi1 -> j * (pi/2), {j, 0, 1}]/2];
{Row[{"Zero quatum coherence signal = ", zq$]},
 Row[{"Double quantum coherence signal = ", dq$}]}
```

```
Out[69]= {Zero quatum coherence signal =
Cos[a13 t1] (Iz1 Cos[a12 t1] + (Mult[Ix1, Iy2] - Mult[Iy1, Ix2]) Sin[a12 t1]) +
(Cos[a12 t1] (Mult[Ix1, Iy3] - Mult[Iy1, Ix3]) -
2 (Mult[Iz1, Ix2, Ix3] + Mult[Iz1, Iy2, Iy3]) Sin[a12 t1]) Sin[a13 t1]

Double quantum coherence signal =
Cos[a13 t1] (Mult[Ix1, Iy2] + Mult[Iy1, Ix2]) Sin[a12 t1] +
(Cos[a12 t1] (Mult[Ix1, Iy3] + Mult[Iy1, Ix3]) +
2 (Mult[Iz1, Ix2, Ix3] - Mult[Iz1, Iy2, Iy3]) Sin[a12 t1]) Sin[a13 t1]}
```

Let's plot the signals by assigning $a_{12} \rightarrow 1$ and $a_{13} \rightarrow 2$

In[105]:=

```

(* Assign values to the dipolar constants *)
dip$ = {a12 → 1., a13 → 2.}; (* in MHz *)
zq$N = zq$ /. dip$;
dq$N = dq$ /. dip$;


(* Plot ZQC *)
modPlotStyle = {Frame → True, Axes → False, FrameStyle → Directive[Thick, Black],
  LabelStyle → Directive[16], FrameLabel → {Style["Time ( $\mu$ s)", 20],
    Style[Row[{"Zero Quantum Coherence\n", "Signal (arbitrary unit)"}], 20}},
  ImageSize → 500, PlotStyle → Directive[Thickness[0.01], Blue],
  FrameTicksStyle → Directive[18], ImagePadding → {{140, 20}, {60, 20}}};
zqps$ = Plot[zq$N /. {Ix1 → 1, Iy1 → 1, Iz1 → 1, Ix2 → 1, Ix3 → 1, Iy2 → 1, Iy3 → 1},
  {t1, 0, 20}, PlotRange → All,
  AxesOrigin → {0, -3},
  FrameTicks → {{{-2, 0, 2}, None}, {{0, 10, 20}, None}}, Evaluate@modPlotStyle];

(* Plot DQC *)
modPlotStyle = {Frame → True, Axes → False, FrameStyle → Directive[Thick, Black],
  LabelStyle → Directive[16], FrameLabel → {Style["Time ( $\mu$ s)", 20],
    Style[Row[{"Double Quantum Coherence\n", "Signal (arbitrary unit)"}], 20}},
  ImageSize → 500, PlotStyle → Directive[Thickness[0.01], Red],
  FrameTicksStyle → Directive[18], ImagePadding → {{140, 20}, {60, 20}}};
dqps$ = Plot[dq$N /. {Ix1 → 1, Iy1 → 1, Iz1 → 1, Ix2 → 1, Ix3 → 1, Iy2 → 1, Iy3 → 1},
  {t1, 0, 20}, PlotRange → All,
  AxesOrigin → {0, -3},
  FrameTicks → {{{-2, 0, 2}, None}, {{0, 10, 20}, None}}, Evaluate@modPlotStyle];

(* Combine the plots *)
plots$ = {{zqps$}, {dqps$}};
Grid[plots$]

```

Out[113]=



The dipolar coupling constant is given by $a_{ij} = \omega_{ij} \times (1 - 3 \cdot \text{Cos}^2[\theta_{ij}])$ and in the solid-state, to obtain a powder-averaged signal, one needs to average the signal over all possible values of θ_{ij} . We set $\omega_{ij} = 52.04 / r_{ij}^3$, where r_{ij} is the distance between spin i and spin j and note that $d\text{Cos}[\theta] = -\text{Sin}[\theta] d\theta$. We further replace the integral with sum over discrete values of θ_{ij} to obtain the powder averaged zero and double quantum coherence signals as follows

```
In[79]:= Clear[zq$PA]
num$ = 40; t1$max = 10.; t1$step = 0.05;
norm$ = Sum[-Sin[θ12] * Sin[θ13] * Sin[θ23],
  {θ12, 0, π, π / num$}, {θ13, 0, π, π / num$}, {θ23, 0, π, π / num$}];
AbsoluteTiming[
  zq$Sum = Sum[-zq$ * Sin[θ12] * Sin[θ13] * Sin[θ23] /. {Ix1 → 1, Iy1 → 1, Iz1 → 1,
    Ix2 → 1, Ix3 → 1, Iy2 → 1, Iy3 → 1} /. {a12 →  $\frac{52.04}{r12^3} * (1 - 3 * \text{Cos}[\theta12]^2)$ ,
    a13 →  $\frac{52.04}{r13^3} * (1 - 3 * \text{Cos}[\theta13]^2)$ , a23 →  $\frac{52.04}{r23^3} * (1 - 3 * \text{Cos}[\theta23]^2)$ } /.
    {r12 → 4., r13 → 3., r23 → 3.}, {θ12, 0, π, π / num$},
    {θ13, 0, π, π / num$}, {θ23, 0, π, π / num$}];
  zq$PA = Table[zq$Sum / norm$, {t1, 0, t1$max, t1$step}];] [[1]]
```

```
Out[82]=
26.5575
```

```
In[83]:= Clear[dq$PA]
num$ = 40;
AbsoluteTiming[
  dq$Sum = Sum[-dq$ * Sin[θ12] * Sin[θ13] * Sin[θ23] /. {Ix1 → 1, Iy1 → 1, Iz1 → 1,
    Ix2 → 1, Ix3 → 1, Iy2 → 1, Iy3 → 1} /. {a12 →  $\frac{52.04}{r12^3} * (1 - 3 * \text{Cos}[\theta12]^2)$ ,
    a13 →  $\frac{52.04}{r13^3} * (1 - 3 * \text{Cos}[\theta13]^2)$ , a23 →  $\frac{52.04}{r23^3} * (1 - 3 * \text{Cos}[\theta23]^2)$ } /.
    {r12 → 4., r13 → 3., r23 → 3.}, {θ12, 0, π, π / num$},
    {θ13, 0, π, π / num$}, {θ23, 0, π, π / num$}];
  dq$PA = Table[dq$Sum / norm$, {t1, 0, t1$max, t1$step}];] [[1]]
```

```
Out[85]=
25.7232
```

Let's plot the powder averaged signals


In[135]:=

```

modPlotStyle = {Frame → True, Axes → False, FrameStyle → Directive[Thick, Black],
  LabelStyle → Directive[16], FrameLabel → {Style["Time ( $\mu$ s)", 24],
    Style[Row[{"Multiple Quantum Coherence\n", "Signal (arbitrary unit)"}], 24]},
  ImageSize → 600, FrameTicksStyle → Directive[18], ImagePadding → {{140, 20}, {60, 20}}};
time$Tab = Table[t1, {t1, 0, t1$max, t1$step}];
ListPlot[{Transpose[{time$Tab, zq$PA}], Transpose[{time$Tab, dq$PA}]], Joined → True,
  PlotRange → All, AxesLabel → {"Time ( $\mu$ s)", "Signal Amplitude"}, AxesOrigin → {0, -0.35},
  FrameTicks → {{-0.25, 0, 0.25, 0.5, 0.75, 1}, None}, {{0, 2, 4, 6, 8, 10}, None}},
  Evaluate@modPlotStyle, PlotLegends → Placed[
    LineLegend[{"Zero Quantum Coherence", "Double Quantum Coherence"}], {0.7, 0.8}],
  PlotStyle → {Directive[Thickness[0.01], Blue], Directive[Thickness[0.01], Red]}]

```

Out[137]:=



To select odd order quantum coherence, we set $\phi_1 = \phi_2$; $\phi_3 = \phi_1 + \pi/2$

```
In[89]:=  $\rho_{\text{OddQC}} = \text{col\$op}[\rho_{5F} /. \{\phi_2 \rightarrow \phi_1, \phi_3 \rightarrow \phi_1 + \pi/2\}, 2, 1]$ 
```

```
Out[89]=
```

```
Iy1 Cos[a12 t1] Cos[a13 t1] Cos[φ1] + 2 Cos[a13 t1] Cos[φ1] Mult[Iz1, Ix2] Sin[a12 t1] +
  2 Cos[a12 t1] Cos[φ1] Mult[Iz1, Ix3] Sin[a13 t1] -
  4 Cos[φ1]3 Mult[Iy1, Ix2, Ix3] Sin[a12 t1] Sin[a13 t1] -
  Ix1 Cos[a12 t1] Cos[a13 t1] Sin[φ1] + 2 Cos[a13 t1] Mult[Iz1, Iy2] Sin[a12 t1] Sin[φ1] +
  2 Cos[a12 t1] Mult[Iz1, Iy3] Sin[a13 t1] Sin[φ1] +
  4 Cos[φ1]2 Mult[Ix1, Ix2, Ix3] Sin[a12 t1] Sin[a13 t1] Sin[φ1] -
  4 Cos[φ1]2 Mult[Iy1, Ix2, Iy3] Sin[a12 t1] Sin[a13 t1] Sin[φ1] -
  4 Cos[φ1]2 Mult[Iy1, Iy2, Ix3] Sin[a12 t1] Sin[a13 t1] Sin[φ1] +
  4 Cos[φ1] Mult[Ix1, Ix2, Iy3] Sin[a12 t1] Sin[a13 t1] Sin[φ1]2 +
  4 Cos[φ1] Mult[Ix1, Iy2, Ix3] Sin[a12 t1] Sin[a13 t1] Sin[φ1]2 -
  4 Cos[φ1] Mult[Iy1, Iy2, Iy3] Sin[a12 t1] Sin[a13 t1] Sin[φ1]2 +
  4 Mult[Ix1, Iy2, Iy3] Sin[a12 t1] Sin[a13 t1] Sin[φ1]3
```

We apply phase cycle to ρ_{OddC} to separate out ± 3 coherence orders from ± 1

```
In[90]:=  $\text{tq\$} = \text{FullSimplify}\left[\text{Sum}\left[(-1)^j * \rho_{\text{OddQC}} /. \phi_1 \rightarrow j * \frac{\pi}{3} + \frac{\pi}{6}, \{j, 0, 5\}\right] / 6\right];$ 
```

```
Row[{"Triple quantum coherence signal = ", tq$}]
```

```
Out[91]=
```

```
Triple quantum coherence signal =
  (Mult[Ix1, Ix2, Ix3] - Mult[Ix1, Iy2, Iy3] - Mult[Iy1, Ix2, Iy3] - Mult[Iy1, Iy2, Ix3])
  Sin[a12 t1] Sin[a13 t1]
```

```
In[92]:= Clear[tq$PA]
```

```
num$ = 40; t1$max = 10.; t1$step = 0.05;
```

```
norm$ = Sum[-Sin[θ12] * Sin[θ13] * Sin[θ23],
  {θ12, 0, π, π/num$}, {θ13, 0, π, π/num$}, {θ23, 0, π, π/num$}];
```

```
AbsoluteTiming[
  tq$Sum = Sum[-tq$ * Sin[θ12] * Sin[θ13] * Sin[θ23] /. {Ix1 → 1, Iy1 → 1, Iz1 → 1,
    Ix2 → 1, Ix3 → 1, Iy2 → 1, Iy3 → 1} /. {a12 →  $\frac{52.04}{r_{12}^3} * (1 - 3 * \text{Cos}[\theta_{12}]^2)$ ,
    a13 →  $\frac{52.04}{r_{13}^3} * (1 - 3 * \text{Cos}[\theta_{13}]^2)$ , a23 →  $\frac{52.04}{r_{23}^3} * (1 - 3 * \text{Cos}[\theta_{23}]^2)$ } /.
    {r12 → 4., r13 → 3., r23 → 3.}, {θ12, 0, π, π/num$},
    {θ13, 0, π, π/num$}, {θ23, 0, π, π/num$}];
  tq$PA = Table[tq$Sum / norm$, {t1, 0, t1$max, t1$step}];][[1]]
```

```
Out[95]=
```

```
20.9583
```


Let's plot the powder averaged triple quantum coherence signal

```

In[96]:= modPlotStyle = {Frame → True, Axes → False, FrameStyle → Black,
  LabelStyle → Directive[16], FrameLabel → {Style["Time ( $\mu$ s)", 24],
    Style[Row[{"Triple Quantum Coherence\n", "Signal (arbitrary unit)"}], 24]},
  ImageSize → 600, PlotStyle → Directive[Thickness[0.01], Orange],
  FrameTicksStyle → Directive[18], ImagePadding → {{140, 20}, {60, 20}}};
time$Tab = Table[t1, {t1, 0, t1$max, t1$step}];
ListPlot[Transpose[{time$Tab, tq$PA}], Joined → True,
  PlotRange → {Automatic, {-0.12, 0.22}}, AxesOrigin → {0, -0.1}, Evaluate@modPlotStyle,
  FrameTicks → {{-0.1, 0, 0.1, 0.2}, None}, {{0, 2, 4, 6, 8, 10}, None}]]

```

Out[98]=



Note that simulating the triple quantum coherence signal using conventional numerical methods would require evaluating matrix exponentials of dimensions 8×8 at $6 \times 201 \times 41 \times 41 = 2,027,286$ points. This number corresponds to density operator evolution over three pulses and three free evolution periods, the time domain sampling, and averaging over θ_{12} and θ_{13} . As the number of spins and/or additional parameters increases, the total number of evaluation points grows rapidly, rendering the computation prohibitively expensive.

It is therefore advantageous to derive an analytical expression for the signal using UniDyn and subsequently export this expression to an efficient numerical platform, such as Python, for performing the required numerical simulations. [1-3]

References

1. Sinha Roy, A., Marohn, J.A. and Freed, J.H., An analysis of double-quantum coherence ESR in an N-spin system: Analytical expressions and predictions. *J. Chem. Phys.* (2024)
 2. Sinha Roy, A., Assafa, T.E., Dzikovski, B., Joshi, N. and Freed, J.H., Rapid Analysis of DEER Signals Including Short Distances. *J. Phys. Chem. Lett.* (2024)
 3. Sinha Roy, A., Tsay, K., Borbat, P.P., Destefano, A., Han, S., Srivastava, M. and Freed, J.H., Detection of Mutation-Induced Conformational Changes in an Intrinsically Disordered Protein by Double Quantum Coherence Electron Spin Resonance Methodology. *J. Am. Chem. Soc.* (2026)
-